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Single use reels and pouches for sterilization are made of:
 - Medical grade paper 60 g/m²
 - Medical laminated film PET12 + PP40 (Polyethylene Terephthalate, 12 µm + Polypropylene, 40 µm)

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SecurEline® pouches and reels are intended to contain medical devices, resist the sterilization process allowing the sterilization of the contents, allow the preservation of the contents in sterile conditions.

SecurEline® pouches and reels are compatible with the following sterilization treatments:

- STEAM
- EtO (Ethylene Oxide)

Pouches and reels are equipped with Type 1 process indicators (EN ISO 11140-1) for Ethylene Oxide and Steam, printed with water-based ink. These provide an indication for the user as they identify whether the barrier system, together with its contents, have already undergone one of the sterilization treatments indicated above. Indicators follow the criterion:

- Steam sterilization = STEAM indicator from pink to brown
- Sterilization with ethylene oxide = EO indicator from blue to yellow / brown

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The sealing on the 4th side for pouches, and on the 3rd and 4th side for reels, can be carried out with all thermosealing machines suitable for the construction materials of SecurEline® products. It is advisable to refer to the International Standard EN 868-5 for the correct setting parameters (temperature, pressure, speed) in order to guarantee the desired minimum resistance according to the sterilization method employed.

The melting temperature of polypropylene is 175°C.

CARATTERISTICHE FISICHE / PHYSICAL CHARACTERISTICS
CARTA MEDICALE 60 g/m² / MEDICAL GRADE PAPER 60 g/m²

PROPERTIES	UNITS	STANDARDS OR METHODS	TYPICAL (TOLERANCE)
SUBSTANCE	g/m ²	ISO 536	60 (+/- 3)
BENDTSEN POROSITY	ml/min	ISO 5636-3	1000 (+/- 250)
AIR PERFORMANCE	µm (Pa.s)	ISO 5636-3	11.4 (+/- 2.8)
BENDTSEN ROUGHNESS FS	ml/min	ISO 8791-2	375 (+/- 125)
BENDTSEN ROUGHNESS WS	ml/min	ISO 8791-2	375 (+/- 125)
PORE SIZE	µm	EN 868-2:2017 (app. D)	21 (+ 14)
THICKNESS	µm	ISO 534	83
TENSILE STRENGTH MD	kN/m	ISO 1924-2	6.4 (- 1.8)
TENSILE STRENGTH CD	kN/m	ISO 1924-2	3.4 (- 1.1)

WET TENSILE STRENGTH MD	kN/m	ISO 3781	2.1 (- 1.1)
WET TENSILE STRENGTH CD	kN/m	ISO 3781	1.1 (- 0.6)
BURST STRENGTH	kPa	ISO 2758	350 (- 120)
WET BURST	kPa	ISO 3689	150 (- 80)
TEARING STRENGTH MD	mN	ISO 1974	600 (- 50)
TEARING STRENGTH CD	mN	ISO 1974	650 (- 100)
COBB TEST (60s)	g/m ²	ISO 535	15 (+ 5)
WATER REPELLENCY	S	EN 868-2:2017 (app. C)	35 (- 15)
FLUORESCENCE	Pts/dm ²	EN 868-2:2017 §4.2.1.7	0

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The product must be stored in its original packaging, at the recommended conditions of temperature (5°C-30°C) and relative humidity (<50%). Avoid exposure to direct heat sources.

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The performance of the product is guaranteed for 5 years from the date of production, if the storage conditions indicated above have been followed.

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Do not reuse. Do not use if package is damaged.

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The material meets the requirements of ISO 10993 – Biological Evaluation of Medical Devices – Part 5: Tests for in vitro Cytotoxicity.

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The material is not subject to 100% inspection. The duty to ensure that the final pack complies with ISO 11607-1 rests with the user of the packaging materials, as they are responsible for validating and controlling requirements such as compatibility with the product, and the conditions necessary to ensure that sterility and performance are maintained during the packaging, sterilization, storage and distribution phases.

EN 868-7:2017 "Adhesive coated paper for low temperature sterilization processes - Requirements and test methods".

EN 868-6:2017 "Paper for low temperature sterilization processes - Requirements and test methods".

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Eline confirms that the product does not have CMR 1A/1B substances in its composition.

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This product is not intentionally formulated with natural rubber latex, including dry natural rubber.

This product is not handled with natural rubber latex gloves.

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Eline products are not intentionally made with phthalates, including but not limited to

(a) Dibutyl phthalate (DBP), benzyl butyl phthalate (BBP), di(2-ethylhexyl)phthalate (DEHP), di-n-octyl phthalate (DnOP), di-isononyl phthalate (DINP), and diisododecyl phthalate (DIDP).

(b) Those listed in the US Consumer Product Safety Improvement Act (CPSIA) of 2008, EU Medical Device Regulation 2017/745 (MDR) referencing Directive 1272/2008/EC Annex VI, RoHS 3, and the European Union Regulation (EC) 1907/2006.

(c) Those listed on California Proposition 65, EU REACH SVHC, and the EPA Phthalate Action Plan lists.

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Concerning Regulation (EC) No. 1907/2006, commonly known as REACH,

(a) Eline products are compliant with its obligations under REACH.

(b) As delivered, Eline products are classified as articles and are exempt from registration requirements.

(c) Eline products do not have reportable substances of very high concern (SVHC).

(d) Eline products comply with Annex XVII restrictions.

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Medical devices must be disposed of in accordance with medical procedures adopted by health care institutions and organizations and local, state or federal laws and regulations.

Sterile Barrier System devices of the SecurEline[®] series are comparable to urban waste and therefore they must be disposed of as urban waste (mixed, non-recyclable).

The PET plastic film has been processed in such a way as to obtain a stratification that makes it non-recyclable. It must therefore be disposed of as non-recyclable waste.

If the operator decouples film from paper, the paper can be recycled in separate collection.

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