

Technical Data Sheet

Manufacturers Name:	Blue Medical B.V.
Manufacturers Address:	Sluispolderweg 79, 1505HJ, Zaandam, The Netherlands
SRN (Single Registration Number):	NL-MF-000010818
Name of the Device (s):	CP Emulating Indicator
Product group code (s):	BM3300206*, BM33002070, BM33002071, BM33002074
Classification:	Non-Medical Device
Manufacturer product category:	CP
Manufacturer product group:	CP Emulating indicator
Market product category:	Chemical indicators
Product color:	Ink turns from blue to pink
Product size:	Indicator strip: 50mm x 20mm
Applicable standards:	ISO-11140-1:2014 (Type 6)
Suitable processes:	STEAM
Stated Values (SV):	BM33002060: 134°C – 3,5 min. / 121°C – 15 min. BM33002061: 134°C – 3,5 min. / 121°C – 15 min. BM33002065: 134°C – 5,3 min. / 121°C – 15 min. BM33002070: 134°C – 7 min. / 121°C – 20 min. BM33002071: 134°C – 7 min. / 121°C – 20 min. BM33002074: 134°C – 18 min.
Instructions for use:	YES, eIFU www.blue-medical.nl/eifu
HS / Customs code:	3822.9000
Shelf life:	4 years
Lead free:	YES
Latex free:	YES

DEHP / CMR free:

YES

Storage conditions:

Humidity limitation 30% - 70% / Temperature limit +10°C - +30°C / Keep dry / Keep away from sunlight

Product symbol markings:

Do not re-use / Do not use if package is damaged / Catalogue number (REF) / Batch code (LOT) / Consult instructions for use / Date of manufacture / Humidity limitation / Keep away from sunlight / Manufacturer / STEAM / Temperature limit / Use-by date / Keep dry

Box content:

See article summary below.

Additional note:

-

Articles valid under this TDS:

Article number	UDI-DI (GTIN)	Description	Box content
BM33002060	8720174715457	CP Emulating Indicator 3.5 min - Type 6 - 500 pcs	500 pcs
BM33002061	8720174715464	CP Emulating Indicator SA 3.5 min - Type 6 - 250 pcs	250 pcs
BM33002065	8720174715471	CP Emulating Indicator 5.3 min - Type 6 - 500 pcs	500 pcs
BM33002070	8720174715488	CP Emulating Indicator 7 min - Type 6 - 500 pcs	500 pcs
BM33002071	8720174715525	CP Emulating Indicator SA 7 min - Type 6 - 250 pcs	250 pcs
BM33002074	8720174715518	CP Emulating Indicator Prion 18 min - Type 6 - 250 pcs	250 pcs

This Technical Data Sheet is issued under the sole responsibility of Blue Medical B.V.. We hereby declare that the product(s) specified above meet the indicated standards. All supporting documentation is retained at the premises of the manufacturer.

Signature:



Place, date (dd.mm.yyyy) of issue and version:

Mr. Alex Welleweerd
Sales Director

Zaandam, 17.02.2023, v8